












OSSIX® PLUS

OSSIX® VOLUMAX

OSSIX™ Bone
BONE GRAFTING MATERIAL

| | Ossifying, Cross-Linked Barrier Membrane | Volumizing, Ossifying Collagen Scaffold | Ossifying Collagen Sponge | |
|---|---|---|---|---|
| Lateral Bone Augmentation | ✓✓ Possible with all bone grafts | ✓✓ (with/without bone graft) Subject to supporting/ non- supporting bone | ✓✓✓ Subject to supporting/ non-supporting bone with OSSIX Plus |  |
| Socket Preservation | ✓✓ with/without bone graft (contained defect) | ✓✓ with/without bone graft | ✓✓✓ no need for membrane |  |
| Crestal Sinus Elevation | n/a | n/a | ✓✓✓ no need for membrane |  |
| Lateral Window Sinus Elevation | ✓✓✓ over lateral window | ✓✓✓ over lateral window | ✓ |  |
| Vertical Augmentation | ✓✓ with bone graft | ✓✓ with bone graft | ✓✓ with membrane |  |
| Preventing / Remediating Buccal Bone Loss | ✓✓ with bone graft | ✓✓✓ without bone graft | ✓✓ with membrane |  |
| 2nd Stage Implant Surgery Dehiscence | ✓✓ with bone graft | ✓✓✓ without bone graft | ✓✓ with membrane |  |
| 2nd Stage Implant Surgery Vertical | ✓✓ with bone graft | ✓✓ without bone graft (Poncho) | ✓✓ with membrane |  |
| Gingival Recession | n/a | ✓ | n/a |  |
| Increase Soft Tissue Volume* | n/a | ✓✓ | n/a |  |
| Mandibular Nerve Repositioning* | n/a | ✓✓ with/without bone graft | n/a |  |